



ALB Antibody, HRP conjugated

Product Code	CSB-PA00180F1Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P85295
Immunogen	Serum Albumin purified from Goat serum (Native Protein)
Raised In	Rabbit
Species Reactivity	Goat
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serum albumin, ALB
Immunogen Species	Goat
Research Area	Others
Target Names	ALB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.