



ALB Antibody, HRP conjugated

Product Code	CSB-PA00260F1Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49822
Immunogen	Serum Albumin purified from Canine serum (Native Protein)
Raised In	Rabbit
Species Reactivity	Canine
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serum albumin (allergen Can f 3), ALB
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Research Area	Cardiovascular
Target Names	ALB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.