



ALB Antibody

Product Code	CSB-PA00260E1Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49822
Immunogen	Serum Albumin purified from Canine serum (Native Protein)
Raised In	Rabbit
Species Reactivity	Canine
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serum albumin (allergen Can f 3), ALB
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Research Area	Cardiovascular
Target Names	ALB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.