




FUS Antibody

Product Code	CSB-PA02704A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35637
Immunogen	Recombinant Human RNA-binding protein FUS protein (1-198AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA-binding protein FUS (75 kDa DNA-pairing protein) (Oncogene FUS) (Oncogene TLS) (POMp75) (Translocated in liposarcoma protein), FUS, TLS
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	FUS

Image



Western blot

All lanes: FUS antibody at 2µg/ml

Lane 1: 293T whole cell lysate


Lane 2: EC109 whole cell lysate

Secondary


Goat polyclonal to rabbit IgG at 1/15000 dilution

Predicted band size: 54 kDa

Observed band size: 60 kDa



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA02704A0Rb at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA02704A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)