



# PITPNM1 Antibody

Product Code	CSB-PA01995A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00562
Immunogen	Recombinant Human Membrane-associated phosphatidylinositol transfer protein 1 protein (1-238AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Membrane-associated phosphatidylinositol transfer protein 1 (Drosophila retinal degeneration B homolog) (Phosphatidylinositol transfer protein, membrane-associated 1) (PITPnm 1) (Pyk2 N-terminal domain-interacting receptor 2) (NIR-2), PITPNM1, DRES9 NIR2 PITPNM
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PITPNM1