



ARF5 Antibody, FITC conjugated

Product Code	CSB-PA00905C0Rb
Abbreviation	ADP-ribosylation factor 5 protein
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P84085
Immunogen	Recombinant Human ADP-ribosylation factor 5 protein (2-180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	GTP-binding protein that functions as an allosteric activator of the cholera toxin catalytic subunit, an ADP-ribosyltransferase. Involved in protein trafficking; may modulate vesicle budding and uncoating within the Golgi apparatus.
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ADP-ribosylation factor 5, ARF5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ARF5