



PSME2 Antibody

Product Code	CSB-PA01134A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UL46
Immunogen	Recombinant Human Proteasome activator complex subunit 2 protein (2-239AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proteasome activator complex subunit 2 (11S regulator complex subunit beta) (REG-beta) (Activator of multicatalytic protease subunit 2) (Proteasome activator 28 subunit beta) (PA28b) (PA28beta), PSME2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSME2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.