




GDI2 Antibody

Product Code	CSB-PA00945A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50395
Immunogen	Recombinant Human Rab GDP dissociation inhibitor beta protein (1-441AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rab GDP dissociation inhibitor beta (Rab GDI beta) (Guanosine diphosphate dissociation inhibitor 2) (GDI-2), GDI2, RABGDIB
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GDI2

Image



Western blot
 All lanes: GDI2 antibody at 2µg/ml
 Lane 1: 293T whole cell lysate
 Lane 2: EC109 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/15000 dilution
 Predicted band size: 51, 46 kDa
 Observed band size: 51 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.