



HMGB1 Antibody, Biotin conjugated

Product Code	CSB-PA01604D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09429
Immunogen	Recombinant Human High mobility group protein B1 protein (8-179AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	High mobility group protein B1 (High mobility group protein 1) (HMG-1), HMGB1, HMG1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HMGB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.