




GM2A Antibody

Product Code	CSB-PA02375A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17900
Immunogen	Recombinant Human Ganglioside GM2 activator protein (32-193AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ganglioside GM2 activator (Cerebroside sulfate activator protein) (GM2-AP) (Sphingolipid activator protein 3) (SAP-3) [Cleaved into: Ganglioside GM2 activator isoform short], GM2A
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	GM2A

Image



Western blot
 All lanes: GM2A antibody at 2µg/ml + Rat kidney tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 21 kDa
 Observed band size: 21 kDa

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.