



# GM2A Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA02375C0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P17900
<b>Immunogen</b>	Recombinant Human Ganglioside GM2 activator protein (32-193AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Ganglioside GM2 activator (Cerebroside sulfate activator protein) (GM2-AP) (Sphingolipid activator protein 3) (SAP-3) [Cleaved into: Ganglioside GM2 activator isoform short], GM2A
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Tags & Cell Markers
<b>Target Names</b>	GM2A
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.