





PFDN5 Antibody

Product Code	CSB-PA01415A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99471
Immunogen	Recombinant Human Prefoldin subunit 5 protein (2-154AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prefoldin subunit 5 (Myc modulator 1) (c-Myc-binding protein Mm-1), PFDN5, MM1 PFD5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PFDN5


Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA01415A0Rb at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA01415A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA01415A0Rb at dilution of 1:100