



PSMA7 Antibody, HRP conjugated

Product Code	CSB-PA01335B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14818
Immunogen	Recombinant Human Proteasome subunit alpha type-7 protein (5-232AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proteasome subunit alpha type-7 (EC 3.4.25.1) (Proteasome subunit RC6-1) (Proteasome subunit XAPC7), PSMA7, HSPC
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	PSMA7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.