




B2M Antibody

Product Code	CSB-PA05786A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61769
Immunogen	Recombinant Human Beta-2-microglobulin protein (22-119AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-2-microglobulin [Cleaved into: Beta-2-microglobulin form pl 5.3], B2M
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	B2M

Image



Western Blot

Positive WB detected in: Mouse thymus tissue

All lanes: B2M antibody at 4µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 14 kDa

Observed band size: 14 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.