



B2M Antibody, HRP conjugated

Product Code	CSB-PA05786B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61769
Immunogen	Recombinant Human Beta-2-microglobulin protein (22-119AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-2-microglobulin [Cleaved into: Beta-2-microglobulin form pl 5.3], B2M
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	B2M
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.