





CST3 Antibody

Product Code	CSB-PA06196A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01034
Immunogen	Recombinant Human Cystatin-C protein (27-146AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:10-1:100
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cystatin-C (Cystatin-3) (Gamma-trace) (Neuroendocrine basic polypeptide) (Post-gamma-globulin), CST3
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	CST3

Image



Western blot
 All lanes: CST3 antibody at 2µg/ml + Human positive serum
 Secondary
 Goat polyclonal to rabbit IgG at 1/15000 dilution
 Predicted band size: 16 kDa
 Observed band size: 16 kDa



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA06196A0Rb at dilution of 1:10

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.