



Recombinant *Chloroflexus aurantiacus* UPF0133 protein Chy400_0558 (Chy400_0558)

Product Code	CSB-EP496906DSW
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B9LJ33
Product Type	Recombinant Protein
Immunogen Species	<i>Chloroflexus aurantiacus</i> (strain ATCC 29364 / DSM 637 / Y-400-fl)
Purity	>85% (SDS-PAGE)
Sequence	MNQRQLMQMA QQMQRQMVKV QEELAATIVE GTAGGGAITV KMNGHREVQS ITISPEVVDP DDVEMLQDLL LVAINDASRK AQLLAEERMQ PLTGGLKGLF
Source	E.coli
Target Names	Chy400_0558
Protein Names	Recommended name: UPF0133 protein Chy400_0558
Expression Region	1-100
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.