



IL1B Antibody, Biotin conjugated

Product Code	CSB-PA06574D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01584
Immunogen	Recombinant Human Interleukin-1 beta protein (117-269AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-1 beta (IL-1 beta) (Catabolin), IL1B, IL1F2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL1B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.