



FKBP1A Antibody, HRP conjugated

Product Code	CSB-PA02325B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62942
Immunogen	Recombinant Human Peptidyl-prolyl cis-trans isomerase FKBP1A protein (2-103AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidyl-prolyl cis-trans isomerase FKBP1A (PPIase FKBP1A) (EC 5.2.1.8) (12 kDa FK506-binding protein) (12 kDa FKBP) (FKBP-12) (Calstabin-1) (FK506-binding protein 1A) (FKBP-1A) (Immunophilin FKBP12) (Rotamase), FKBP1A, FKBP1 FKBP12
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	FKBP1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.