



CCL2 Antibody, HRP conjugated

Product Code	CSB-PA05865B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13500
Immunogen	Recombinant Human C-C motif chemokine 2 protein (24-99AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-C motif chemokine 2 (HC11) (Monocyte chemoattractant protein 1) (Monocyte chemotactic and activating factor) (MCAF) (Monocyte chemotactic protein 1) (MCP-1) (Monocyte secretory protein JE) (Small-inducible cytokine A2), CCL2, MCP1 SCYA2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CCL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.