



PPP1R11 Antibody, FITC conjugated

Product Code	CSB-PA03834C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60927
Immunogen	Recombinant Human E3 ubiquitin-protein ligase PPP1R11 protein (1-126AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase PPP1R11 (EC 2.3.2.27) (Hemochromatosis candidate gene V protein) (HCG V) (Protein phosphatase 1 regulatory subunit 11) (Protein phosphatase inhibitor 3), PPP1R11, HCGV TCTE5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PPP1R11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.