





CSNK2B Antibody

Product Code	CSB-PA01655A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P67870
Immunogen	Recombinant Human Casein kinase II subunit beta protein (1-215AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Casein kinase II subunit beta (CK II beta) (Phosvitin) (Protein G5a), CSNK2B, CK2N G5A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CSNK2B

Image



Immunohistochemistry of paraffin-embedded human ovarian cancer using CSB-PA01655A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA01655A0Rb at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.