



PMP2 Antibody, FITC conjugated

Product Code	CSB-PA07854C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02689
Immunogen	Recombinant Human Myelin P2 protein (1-132AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myelin P2 protein (Peripheral myelin protein 2), PMP2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PMP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.