




HADHB Antibody

Product Code	CSB-PA03035A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55084
Immunogen	Recombinant Human Trifunctional enzyme subunit beta, mitochondrial protein (35-283AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Trifunctional enzyme subunit beta, mitochondrial (TP-beta) [Includes: 3-ketoacyl-CoA thiolase (EC 2.3.1.16) (Acetyl-CoA acyltransferase) (Beta-ketothiolase)], HADHB
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HADHB

Image



Immunohistochemistry of paraffin-embedded human prostate tissue using CSB-PA03035A0Rb at dilution of 1:50