



DAD1 Antibody, Biotin conjugated

Product Code	CSB-PA03124D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61803
Immunogen	Recombinant Human Dolichyl-diphosphooligosaccharide--protein glycosyltransferase subunit DAD1 protein (1-113AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dolichyl-diphosphooligosaccharide--protein glycosyltransferase subunit DAD1 (Oligosaccharyl transferase subunit DAD1) (Defender against cell death 1) (DAD-1), DAD1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DAD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.