



PGK1 Antibody, HRP conjugated

Product Code	CSB-PA00035B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00558
Immunogen	Recombinant Human Phosphoglycerate kinase 1 protein (25-417AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphoglycerate kinase 1 (EC 2.7.2.3) (Cell migration-inducing gene 10 protein) (Primer recognition protein 2) (PRP 2), PGK1, PGKA
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PGK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.