



PGK1 Antibody, FITC conjugated

Product Code	CSB-PA00035C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00558
Immunogen	Recombinant Human Phosphoglycerate kinase 1 protein (25-417AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphoglycerate kinase 1 (EC 2.7.2.3) (Cell migration-inducing gene 10 protein) (Primer recognition protein 2) (PRP 2), PGK1, PGKA
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PGK1