




MP68 Antibody

Product Code	CSB-PA00675A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56378
Immunogen	Recombinant Human 6.8 kDa mitochondrial proteolipid protein (1-58AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MP68

Image



Western Blot

Positive WB detected in: PC-3 whole cell lysate, 293T whole cell lysate, THP-1 whole cell lysate, Jurkat whole cell lysate, A549 whole cell lysate, SH-SY5Y whole cell lysate, HeLa whole cell lysate

All lanes: MP68 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 7, 9 kDa

Observed band size: 7 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.