



IL1A Antibody

Product Code	CSB-PA06565A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01583
Immunogen	Recombinant Human Interleukin-1 alpha protein (113-271AA)
Raised In	Rabbit
Species Reactivity	Human
Specificity	This antibody will recognise the precursor and mature form of IL-1a (i.e 31kDa and 17kDa forms). The antiserum does not recognize human IL-1b or Mouse or Rabbit IL-1a. This antibody is also useful for neutralization of human and primate IL-1a activity in bioassays. It does not neutralize the biological activity IL-1b. It does not neutralize the biological activity of mouse, rat or rabbit IL-1a. It will neutralize primate/human IL-a.
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-1 alpha (IL-1 alpha) (Hematopoietin-1), IL1A, IL1F1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.