



IL1A Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA06565C0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P01583 |
| Immunogen | Recombinant Human Interleukin-1 alpha protein (113-271AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Specificity | This antibody will recognise the precursor and mature form of IL-1a (i.e 31kDa and 17kDa forms). The antiserum does not recognize human IL-1b or Mouse or Rabbit IL-1a. This antibody is also useful for neutralization of human and primate IL-1a activity in bioassays. It does not neutralize the biological activity IL-1b. It does not neutralize the biological activity of mouse, rat or rabbit IL-1a. It will neutralize primate/human IL-a. |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Interleukin-1 alpha (IL-1 alpha) (Hematopoietin-1), IL1A, IL1F1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | IL1A |