



RRAGA Antibody, FITC conjugated

Product Code	CSB-PA03465C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L523
Immunogen	Recombinant Human Ras-related GTP-binding protein A protein (1-247AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-related GTP-binding protein A (Rag A) (RagA) (Adenovirus E3 14.7 kDa-interacting protein 1) (FIP-1), RRAGA
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	RRAGA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.