




# SNRPD2 Antibody

<b>Product Code</b>	CSB-PA02635A0Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P62316
<b>Immunogen</b>	Recombinant Human Small nuclear ribonucleoprotein Sm D2 protein (1-118AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Small nuclear ribonucleoprotein Sm D2 (Sm-D2) (snRNP core protein D2), SNRPD2, SNRPD1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	SNRPD2

## Image



Western blot  
All lanes: SNRPD2 antibody at 2µg/ml + 293T whole cell lysate  
Secondary  
Goat polyclonal to rabbit IgG at 1/15000 dilution  
Predicted band size: 14, 13 kDa  
Observed band size: 14 kDa