



PHGDH Antibody, Biotin conjugated

Product Code	CSB-PA03255D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43175
Immunogen	Recombinant Human D-3-phosphoglycerate dehydrogenase protein (2-251AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	D-3-phosphoglycerate dehydrogenase (3-PGDH) (EC 1.1.1.95) (2-oxoglutarate reductase) (EC 1.1.1.399) (Malate dehydrogenase) (EC 1.1.1.37), PHGDH, PGDH3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PHGDH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.