




ARPC2 Antibody

Product Code	CSB-PA03365A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15144
Immunogen	Recombinant Human Actin-related protein 2/3 complex subunit 2 protein (1-250AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Actin-related protein 2/3 complex subunit 2 (Arp2/3 complex 34 kDa subunit) (p34-ARC), ARPC2, ARC34
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ARPC2

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA03365A0Rb at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.