



IFITM1 Antibody, Biotin conjugated

Product Code	CSB-PA01204D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13164
Immunogen	Recombinant Human Interferon-induced transmembrane protein 1 protein (1-36AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon-induced transmembrane protein 1 (Dispanin subfamily A member 2a) (DSPA2a) (Interferon-induced protein 17) (Interferon-inducible protein 9-27) (Leu-13 antigen) (CD antigen CD225), IFITM1, CD225 IFI17
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IFITM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.