



CXCL5 Antibody, HRP conjugated

Product Code	CSB-PA06155B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42830
Immunogen	Recombinant Human C-X-C motif chemokine 5 protein (37-110AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-X-C motif chemokine 5 (ENA-78(1-78)) (Epithelial-derived neutrophil-activating protein 78) (Neutrophil-activating peptide ENA-78) (Small-inducible cytokine B5) [Cleaved into: ENA-78(8-78); ENA-78(9-78)], CXCL5, ENA78 SCYB5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	CXCL5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.