



ERP29 Antibody, Biotin conjugated

Product Code	CSB-PA05345D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30040
Immunogen	Recombinant Human Endoplasmic reticulum resident protein 29 protein (40-251AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endoplasmic reticulum resident protein 29 (ERp29) (Endoplasmic reticulum resident protein 28) (ERp28) (Endoplasmic reticulum resident protein 31) (ERp31), ERP29, C12orf8 ERP28
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ERP29
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.