




VOPP1 Antibody

Product Code	CSB-PA03545A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96AW1
Immunogen	Recombinant Human Vesicular, overexpressed in cancer, prosurvival protein 1 protein (23-172AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:1000-1:5000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vesicular, overexpressed in cancer, prosurvival protein 1 (EGFR-coamplified and overexpressed protein) (ECop) (Glioblastoma-amplified secreted protein) (Putative NF-kappa-B-activating protein 055N), VOPP1, ECOP GASP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	VOPP1

Image



Western blot

All lanes: VOPP1 antibody at 2µg/ml

Lane 1: EC109 whole cell lysate


Lane 2: 293T whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/15000 dilution

Predicted band size: 20, 19, 18 kDa

Observed band size: 50, 70, 80 kDa



Immunofluorescent analysis of HepG2 cells using CSB-PA03545A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.