



VOPP1 Antibody, Biotin conjugated

Product Code	CSB-PA03545D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96AW1
Immunogen	Recombinant Human Vesicular, overexpressed in cancer, prosurvival protein 1 protein (23-172AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vesicular, overexpressed in cancer, prosurvival protein 1 (EGFR-coamplified and overexpressed protein) (ECop) (Glioblastoma-amplified secreted protein) (Putative NF-kappa-B-activating protein 055N), VOPP1, ECOP GASP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	VOPP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.