




RPL15 Antibody

Product Code	CSB-PA02435A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61313
Immunogen	Recombinant Human 60S ribosomal protein L15 protein (2-204AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RPL15

Image



Western Blot

Positive WB detected in: A549 whole cell lysate(30µg), HeLa whole cell lysate(30µg), MCF7 whole cell lysate(30µg)

All lanes: RPL15 antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/20000 dilution

Predicted band size: 25, 17 kDa

Observed band size: 17 kDa

Exposure time: 60s

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.