



RPL15 Antibody, HRP conjugated

Product Code	CSB-PA02435B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61313
Immunogen	Recombinant Human 60S ribosomal protein L15 protein (2-204AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	60S ribosomal protein L15 (Large ribosomal subunit protein eL15), RPL15, EC45
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RPL15
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.