



Recombinant Human Protachykinin-1 (TAC1)

| | |
|--------------------------|--|
| Product Code | CSB-BP023062HU |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P20366 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MKILVALAVFFLVSTQLFAEEIGANDDLNYWSDWYDSDQIKEELPEPFEHLLQR IARRPKPQQFFGLMGKRDADSSIEKQVALLKALYGHGQISHKRHKTDTSFVGLM GKRALNSVAYERSAMQNYERRR |
| Source | Baculovirus |
| Target Names | TAC1 |
| Protein Names | Recommended name: Protachykinin-1 Alternative name(s): PPT Cleaved into the following 5 chains: 1. Substance P 2. Neurokinin A Short name= 3. NKA Alternative name(s): Neuromedin L Substance K Neuropeptide K Short nam |
| Expression Region | 1-129 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Cytoplasmic domain |
| Target Details | This gene encodes four products of the tachykinin peptide hormone family, substance P and neurokinin A, as well as the related peptides, neuropeptide K and neuropeptide gamma. These hormones are thought to function as neurotransmitters which interact with nerve receptors and smooth muscle cells. They are known to induce behavioral responses and function as vasodilators and secretagogues. Multiple transcript variants encoding different isoforms have been found for this gene. |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |