




RPS18 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA01645A0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P62269 |
| Immunogen | Recombinant Human 40S ribosomal protein S18 protein (2-152AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, WB; Recommended dilution: WB:1:1000-1:5000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | 40S ribosomal protein S18 (Ke-3) (Ke3) (Small ribosomal subunit protein uS13), RPS18, D6S218E |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling |
| Target Names | RPS18 |

Image



Western blot
 All lanes: RPS18 antibody at 2µg/ml + 293T whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/15000 dilution
 Predicted band size: 18 kDa
 Observed band size: 75 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.