



TIMP2 Antibody, Biotin conjugated

Product Code	CSB-PA07405D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16035
Immunogen	Recombinant Human Metalloproteinase inhibitor 2 protein (30-220AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Metalloproteinase inhibitor 2 (CSC-21K) (Tissue inhibitor of metalloproteinases 2) (TIMP-2), TIMP2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	TIMP2