



VEGFC Antibody, Biotin conjugated

Product Code	CSB-PA07545D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49767
Immunogen	Recombinant Human Vascular endothelial growth factor C protein (112-227AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vascular endothelial growth factor C (VEGF-C) (Flt4 ligand) (Flt4-L) (Vascular endothelial growth factor-related protein) (VRP), VEGFC
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	VEGFC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.