




TAX1BP3 Antibody

Product Code	CSB-PA04445A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14907
Immunogen	Recombinant Human Tax1-binding protein 3 protein (1-124AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tax1-binding protein 3 (Glutaminase-interacting protein 3) (Tax interaction protein 1) (TIP-1) (Tax-interacting protein 1), TAX1BP3, TIP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TAX1BP3

Image



Western blot

All lanes: TAX1BP3 antibody at 2µg/ml

Lane 1: EC109 whole cell lysate

Lane 2: 293T whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/15000 dilution

Predicted band size: 14 kDa

Observed band size: 70 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.