



SYNCRIP Antibody, FITC conjugated

Product Code	CSB-PA03334C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60506
Immunogen	Recombinant Human Heterogeneous nuclear ribonucleoprotein Q protein (2-191AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heterogeneous nuclear ribonucleoprotein Q (hnRNP Q) (Glycine- and tyrosine-rich RNA-binding protein) (GRY-RBP) (NS1-associated protein 1) (Synaptotagmin-binding, cytoplasmic RNA-interacting protein), SYNCRIP, HNRPQ NSAP1
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	SYNCRIP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.