



# GABARAP Antibody, Biotin conjugated

Product Code	CSB-PA05484D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95166
Immunogen	Recombinant Human Gamma-aminobutyric acid receptor-associated protein (2-117AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gamma-aminobutyric acid receptor-associated protein (GABA(A) receptor-associated protein) (MM46), GABARAP, FLC3B
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GABARAP