




dgcQ Antibody

Product Code	CSB-PA08524A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P76330
Immunogen	Recombinant Escherichia coli Probable diguanylate cyclase DgcQ protein (381-564AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable diguanylate cyclase DgcQ (DGC) (EC 2.7.7.65) (Cellulose synthesis regulatory protein), dgcQ, yedQ
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	dgcQ

Image



Western Blot

Positive WB detected in: Rosseta bacteria lysate at 80µg, 40µg, 20µg, 10µg

All lanes: yedQ antibody at 2.5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

predicted band size: 65 kDa

observed band size: 65 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.