





AP2S1 Antibody

Product Code	CSB-PA02604A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53680
Immunogen	Recombinant Human AP-2 complex subunit sigma protein (1-142AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	AP-2 complex subunit sigma (Adaptor protein complex AP-2 subunit sigma) (Adaptor-related protein complex 2 subunit sigma) (Clathrin assembly protein 2 sigma small chain) (Clathrin coat assembly protein AP17) (Clathrin coat-associated protein AP17) (HA2 17 kDa subunit) (Plasma membrane adaptor AP-2 17 kDa protein) (Sigma2-adaptin), AP2S1, AP17 CLAPS2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AP2S1

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA02604A0Rb at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA02604A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.